



Slide storage boxes

Manufactured in polystyrene, these economical boxes are ideal for storing or transport 25, 50 or 100 76 x 26 mm slides.

50 and 100 boxes include numbered sheets with a labelling in both the bottom of the case and in the lid (25 box does not include lid label).

They have a hinged lid for safe closing and are easy opening. Available in beige.

Automatic slide dispenser

Made of ABS and PS, this dispenser can hold up to 50 standard slides. Slides are dispensed individually by rotating the knurled knob on either side of the moulded base.



Microscopes Slides

Good quality, low cost glass slides (50 pack) for routine lab use.

Coverslips

To compliment the microscopes slides we also supply coverslips Available in boxes of 100.

See price list for sizes stocked.

